




Kineticopro®
Total Water Care

Kineticopro Water Test Kit Reference Guide

1. You will receive a box labeled 'Water Sample Test Kit' from Kineticopro. This will include a return shipment label once you have completed the instructions inside.
2. Inside the box you will find 2 bottles and a 5-Way Water Check kit. Please remove items from the box.
3. The 1st bottle, labeled 'Unfiltered Water', will be used to capture a sample of water from water source that IS NOT currently being treated. (ex: bathroom sink, kitchen sink, etc.).
4. The 2nd bottle, labeled 'Filtered Water', will be used to capture a water sample from a treated water source, if applicable (ex: water from fountain beverage machine, filtered drink water station, etc.).
5. Please utilize the 5-Way Water Check kit to identify on-site levels of chlorine. Simply follow the instructions on the back of the packet and be sure to mark the chlorine levels for each sample on the bottle labels.
6. When filling up the water bottles please be sure to follow the instructions provided on the side of each water bottle, which can also be found inside the box.
7. After water bottles are filled up per the instructions, please ensure lids are screwed on tightly and put back into the box provided. Then adhere return label on the outside of the box and call UPS to make pick-up arrangements OR drop off at your nearest UPS store.

All bottles need to be returned to the following address:

Kineticopro Incorporated
3980 Lakefield Ct.,
Suwanee GA, 30024



1.800.321.5022
www.Kineticopro.com
WaterCare@Kineticopro.com